

The Importance of Heptahelical Receptors in "Drugable Space"

"the most fruitful basis for the discovery of a new drug is to start with an old drug" – Sir James W. Black

Only about 6% of NMEs (new molecular entities) target a previously undrugged target

A numbers game:

- 21,000 drug products
- 1,357 unique drugs
- 1,204 "small molecules"
- 803 oral
- 421 parenteral
- 275 topical
- 166 biologicals

For approved drugs:

- 324 molecular targets
- 266 human targets
- 207 human targets for "small molecules"

Genome perspective:

- 1,620 out of 22,218 human genes have explicit disease associations – 105 are drug targets

Rates of Innovation (past 20 years):

- First-against-target – 5.3/year
- New target families – 1.9/year

Rhodopsin-like GPCRs

27%

726 heptahelical receptors encoded by the human genome. ~200 are "orphan" receptors, most expressed in CNS

Overington, J.P. et al., Nat Revs Drug Dis 5:993 (2006), Insight Lecture at Rockefeller University, P.N. Goodfellow, Oct. 22, 2007

Why Is the Study of Heptahelical Receptors So Challenging?

- No modular structure**
We can't study individual domains one at a time
- Membranes are required**
- No orthologues in lower organisms:**
Only eukaryotic organisms have heptahelicals
Bacteriorhodopsins don't count
Can't express in *E. coli*
- Formalisms used to describe GPCR function are 50 years out of date (i.e., classical pharmacology)**
- Low abundance in Nature:**
It can be argued that only 3 GPCRs have been purified in functional form from tissue

Trends in Pharmacological Sciences

Escaping the flatlands: new approaches to study the dynamics of GPCR signaling complexes

See Huber & Sakmar, TIPS, July 2011

Structural Biology of the GPCR Signaling Complex: A Series of Crystallographic Snapshots

1. 1993
2. 1994
3. 1996
4. 2000
5. 2007
6. 2008
7. 2011
8. 2011

β 2AR-Gs complex – Rasmussen *et al.*, *Nature* (2011)
Schwartz & Sakmar *Nature* (2011)

But wait . . . the dynamics of the GPCR-G protein cycle is extremely complicated

- Crystallography has its limitations -- static portraits, often in contrived poses
- Our aim is to complement high-resolution structures with biochemical and biophysical data obtained by other approaches

How will we do that? New technologies are needed . . .

Huber & Sakmar *TIPS* (2011)

New Technologies to Study Receptor Signaling

PART 1
Genetic Encoding of Unnatural Amino Acids in Expressed GPCRs Using Amber Codon Suppression Technology
“Site-directed Unnatural Amino (UAA) Acid Mutagenesis”
(with Shixin Ye, Thomas Huber, Saranga Naganathan, Amy Grunbeck, Sarmistha Ray-Saha, Kelly Daggett)

PART 2
Applications for UAA Mutagenesis and Converging Novel Technologies To Study Structural Dynamics of Receptor “Signalosomes”
FTIR spectroscopy of azido groups, site-directed photocrosslinking, bioorthogonal labeling, oriented-tethered bilayers & SMD-TIRF
(Adam Knepp, Thomas Huber, Saranga Naganathan, Amy Grunbeck, Kelly Daggett, Sarmistha Ray-Saha, Louise Valentin-Hansen)

Expanding the Genetic Code in Mammalian Cells Using Amber Codon Suppression

Engineered E. coli aminoacyl-tRNA synthetase
orthogonal suppressor tRNA (tRNA^{Yam})
Mammalian Cell

UGA: A Third Nonsense Triple in the Genetic Code

mention Alan Garen

Suppressor tRNA in Various Plasmids

5'AGCGCTCCGGT777TGTGTGCTGAACCT3'
3'CAGGGGACGCCGACACAGTACAGGT5'

Nucleic Acids Research, 2002, Vol. 30, 4692-4699
PNAS, 2004, Vol. 101, 8883-8887
Molecular Cell Biology, 1996, Vol. 16, 907-913

Suppressor tRNA derived from *Bacillus stearothermophilus* tyr-tRNA

with Caroline Koehler & Tom RajBhandary, M.I.T.

Using a luciferase reporter system to validate and optimize the mutagenesis scheme

Evaluate tRNA constructs

A.

wt-Fluc	-	p0blunt.3XYam	p0blunt.Yam	pSVB.fMam	pZeo.Yam	pSVB.Yam
---------	---	---------------	-------------	-----------	----------	----------




B.

tRNA construct	Bioluminescence (photons/sec)
no suppressor tRNA	~0
p0blunt.3Yam	~1
p0blunt.Yam	~2
pSVB.fMam	~5
pZeo.Yam	~10
pSVB.Yam	~25

C.

tRNA	F/ug	F%
p0blunt.3XYam	2.87X10 ³	0.01%
pZeo.Yam	4.84X10 ⁵	14.03%
pSVB.Yam	10.11X10 ⁶	29.34%

3.5ug Fluc/3.5ug tRNA/0.35ug RS
10cm plate co-transfection



PART 2 – Applications for Genetic Encoding of Non-natural Amino Acids in Expressed GPCRs (Unnatural Amino Acid Mutagenesis)*

- Tracking receptor activation using FTIR spectroscopy
- Unnatural amino acids as targeted photocrosslinkers
- Bioorthogonal labeling chemistries: e.g. Staudinger Ligation
 - Antibody Epitope Tagging in the native bilayer environment
 - Stoichiometric fluorescence labeling for single-molecule detection
- Single-molecule detection TIRF imaging

Tyrosine analogues amenable to genetic encoding

• chemical handle for labeling

• chemical handle for labeling



• photocrosslinker


• chemical handle for labeling

• photocrosslinker

• infrared probe

*Including Rhodopsin & CXCR4/CCR5 as model systems









PART 2 – Applications for Genetic Encoding of Non-natural Amino Acids in Expressed GPCRs (Unnatural Amino Acid Mutagenesis)*







- Tracking receptor activation using FTIR spectroscopy
- Unnatural amino acids as targeted photocrosslinkers
- Bioorthogonal labeling chemistries: e.g. Staudinger Ligation
Antibody Epitope Tagging in the native bilayer environment
Stoichiometric fluorescence labeling for single-molecule detection
- Single-molecule detection TIRF imaging

Tyrosine analogues amenable to genetic encoding

• chemical handle for labeling	• chemical handle for labeling	• chemical handle for labeling	• chemical handle for labeling
	• photocrosslinker	• photocrosslinker	• photocrosslinker
		• infrared probe	

*Including Rhodopsin & CXCR4/CCR5 as model systems





PART 2 – Applications for Genetic Encoding of Non-natural Amino Acids in Expressed GPCRs (Unnatural Amino Acid Mutagenesis)*

- Tracking receptor activation using FTIR spectroscopy
- Unnatural amino acids as targeted photocrosslinkers
- Bioorthogonal labeling chemistries: e.g. Staudinger Ligation
Antibody Epitope Tagging in the native bilayer environment
Stoichiometric fluorescence labeling for single-molecule detection
- Single-molecule detection TIRF imaging

Tyrosine analogues amenable to genetic encoding

p-acetyl-L-phenylalanine (AcP)	p-benzoyl-L-phenylalanine (BzP)	p-azido-L-phenylalanine (AzP)	p-O-propargyl-L-phenylalanine (PpP)
• chemical handle for labeling	• chemical handle for labeling	• chemical handle for labeling	• chemical handle for labeling
	• photocrosslinker	• photocrosslinker	• infrared probe

*Including Rhodopsin & CXCR4/CCR5 as model systems

Single-molecule detection (SMD) using TIRF (total internal reflectance fluorescence) microscopy and microfluidics

Self assembly of an immobilized membrane bilayer containing oriented recombinant expressed CCR5

- CCR5 c-term 104 mAb epitope
- POPC bilayer
- biotinylated PEG 2000 DSPE
- biotinylated 104 mAb
- streptavidin
- biotinylated BSA
- C-18 silane glass

Chemokines (e.g. SDF1α) can be FL labeled and flowed through TIRF cell.

Toward SMD Fluorescence Microscopy of GPCRs in Live Cell Membranes

A stable oriented "tethered" bilayer can be assembled in a flow-cell "chip" for SMD-TIRF microscopy

Huber & Sakmar (2011) TIPS

Back to the Vertebrate Visual System . . .






Some thoughts on the supramolecular structure of the rod outer segment disc membrane



"How does the G protein (transducin) find an active receptor (R)? And then what?"*

(Adam Knepp, Thomas Huber – With Xavier Periole and Siewert-Jan Marrink, Groningen, The Netherlands)

Model of a heptahelical receptor-G protein complex ("Signalosome") in a membrane bilayer

But How Does the G Protein Find an Active Receptor in the Bilayer?





The Rockefeller University www.SakmarLab.org

2011

Head of Laboratory: Thomas F. Sakmar, M.D.
Robert M. Rosenfeld Professor
Box 25, 1230 York Avenue
New York, NY 10021-6300

817-527-8680 (M-F, 9AM-5PM)
212-327-8686 Office
212-327-1964 Fax
212-327-8236 Cell Laboratory
212-327-8883 Biomechanics Lab

Lab members (L-R): robertros@rockefeller.edu
Tom Palmer, chakrapan
Jennifer Priddy, jpriddy@rockefeller.edu
Nehal Shah, nshah@rockefeller.edu
De Tian, dtian@rockefeller.edu
Scott Lippell, slippell@rockefeller.edu
Levin Sabatini-Duncan, dsabatini@rockefeller.edu
Kaitlyn Kagan, kagank@rockefeller.edu
Thomas Sakmar, tsakmar@rockefeller.edu
Katie Hooper, kahooper@rockefeller.edu
Srinivasa Ray Saha, srinivasa@rockefeller.edu
Adam Rosen, arosen@rockefeller.edu
Sangeeta Reganathan, sreganath@rockefeller.edu
Amy Kravitz, akravitz@rockefeller.edu
Thomas Heide, ttheide@rockefeller.edu
Wendy Wang, wwang@rockefeller.edu
Amy Tang, atang@rockefeller.edu

Check out our new website
SakmarLab.org
and come to see your fave
(don't be shy for that!)

Cover art by Karina Aberg Sakmar

Trends in Pharmacological Sciences

Escaping the flatlands:
new approaches to study the
dynamics of GPCR signaling complexes

Cell